

norma

Diagnostic Solutions

icon

A Smart Hematology Analyzer

- Smallest sample volume makes it ideal for pediatric use
- Smallest reagent consumption
- A4 paper size footprint
- ER mode for use at the point-of-care
- NormaAccess – Remote support tool



www.normadiagnostika.com

TECHNICAL DATA

22 Parameters:

Histograms:

Sample volume:

Sample tube:

Throughput:

Measuring technology:

Storage capacity:

User interface:

Dimensions (HxWxD):

Net Weight:

WBC, LYM, MID, GRA, LYM%, MID%, GRA%,
RBC, HGB, MCV, HCT, MCH, MCHC, RDWcv, RDWsd,
PLT, MPV, PCT, PDWcw, PDWsd, P-LCR%, P-LCC

WBC, RBC, PLT

9.6 μ L whole blood (with EDTA)

closed or open vials

60 tests per hour

volumetric impedance method with microfluidic technology

40,000 test results incl. histograms

10.1 inch, 1280 x 800 resolution LCD with touch screen

28 cm x 22 cm x 32 cm

8.7 kg

Norma Icon is ideal for doctors' offices, emergency departments, small to medium size labs and polyclinics.



- Microfluidic technology
- Patented, innovative solutions
- Trendsetting design



Shortlisted for the Senetics Innovation Award 2014
Category "Most Innovative Product in the Health Care Sector"

reddot award 2014
winner

Norma Instruments Zrt.

Papírgyár u. 58-59., 1038 Budapest
Hungary

Phone: +36 1 815 4370

E-mail: office@normadiagnostika.com



Norma Instruments Zrt. reserves the right without notice to make any technical changes it considers necessary or desirable.

icon