

Norma Autoloader

Automated sample processing for better workflow and time efficiency

- The only retrofittable autosampler on the market
- Suitable for medium and big laboratories
- Built-in sample mixer and barcode reader
- Space-saving compact size
- Sample capacity: 50 vials
- Emergency sample option



norma

Technical specifications

Reference number	NICA-AL50
Compatible with	Icon-5
Accessories	Normaboard, 2 interchangeable trays, power supply
Capacity	50 samples/tray
Supported vial types	Greiner Vacuette, BD Vacutainer, BD MAP
Operation modes	Full sequence operation (full tray scan) – can be paused at any time Random access mode (sample positions can be selected on the analyzer)
Connectivity	2 x USB ports, built-in barcode reader
Dimensions (HxWxD)	340 x 220 x 310 mm
Normaboard Dimensions (HxW)	435 x 375 mm
Weight	8.5 kg
Throughput	60 tests / hour with Icon-5
Power	External power supply, universal, 24 VDC, 2,7 A (100-240 VAC 50-60Hz), non-grounded
Power consumption	Maximum 60W
Other features	▪ Emergency samples can be measured directly on analyzer's closed vial port ▪ All operations can be paused at any time ▪ Connected via USB to the analyzer ▪ Normaboard is used to align analyzer with Autoloader ▪ Built-in sample mixer ▪ Gripper is designed to grab only vials with cap ▪ Tube detection by mechanical force feedback ▪ Removable sample tray with identification



Norma Instruments Zrt.

Papírgyár u. 58-59., 1038 Budapest, Hungary
+36 1 815 4370; office@normadiagnostika.com

Norma Instruments Zrt. reserves the right without notice to make any technical changes it considers necessary or desirable.